

DNA-Flow technology


Molecular diagnostic kits
based on one-step multiplex
PCR and reverse hybridisation
on a membrane.



DNA-Flow technology




Speed



Sensitivity



Multiplex




Flexibility

Vitro has a wide range of molecular diagnostic systems for the simultaneous detection of pathogen panels and markers, based on multiplex PCR or RT/PCR amplification and hybridisation using DNA-FLOW technology.

It is an innovative technology that uses automated flow-through reverse hybridisation on porous membrane-based microarrays that contain oligonucleotides specific to each marker. The target DNA molecules cross the membrane through a vertical vacuum flow, facilitating their interaction with the specific and complementary DNA probes immobilised on the porous membrane in a three-dimensional environment, unlike conventional dot blot technology in which hybridisation occurs passively. This accelerates the process and increases sensitivity, reducing reaction volumes and overall processing time from several hours to a few minutes.


Once binding has occurred between the specific amplicons and their corresponding probes, this is visualised by a colorimetric immunoenzymatic reaction with Streptavidin-Alkaline Phosphatase and the BCIP/NBT chromogen, forming insoluble precipitates on the membrane in those areas where hybridisation has occurred. The process is performed automatically in the hybriSpot™ platform and the results are analysed with the hybriSoft™ software.

Time: Several Hours




Conventional
“Passive”
Hybridization

Performed with Orbital Shaker




Passive diffusion is a slow and random process




- Complimentary DNA
- Non-complimentary DNA
- Immobilized Oligoprobe

Time: 35min



Flow-Through
“Active”
Hybridization



Lamina flow directs target DNA to DNA probe



HybriSpot Systems



REFERENCE: VIT-HS24A



REFERENCE: VIT-HS12A



Automatic systems

- ✓ Fully automated amplification, denaturation and hybridisation.
- ✓ Minimal handling.
- ✓ UV lamp for residual DNA/RNA decontamination.
- ✓ Automatic scanning of reagent barcodes.
- ✓ Barcode scanner compatible with any type of label.
- ✓ HybriSpot software for sample management, automatic image capture and automatic analysis of results.
- ✓ Connects to the LIS in the laboratory.



REFERENCE: VIT-HS12


Manual system

- ✓ Hybridisation process guided via the instrument's display.
- ✓ 12 samples in 45 minutes.
- ✓ HybriSpot software for sample management, automatic image capture, automatic analysis of results.
- ✓ Compact size (350x230x200 mm) and weighs 4 kg.
- ✓ Suitable for laboratories with small numbers of samples.

KITS

HPV DIRECT FLOW CHIP KIT

Detects and identifies HPV virus and simultaneously genotype 35 types of low- and high-risk HPV types



| | | | | | | | | |
|----|----|----|-------|----|----|----|-------|----|
| B | 33 | 58 | 42 | 71 | 16 | 52 | B | |
| B | 35 | 59 | 43 | 72 | 18 | 53 | 6 | 69 |
| C | 39 | 66 | 44/55 | | 26 | 56 | 11 | 70 |
| U | 45 | 68 | 54 | 84 | 31 | 58 | 40 | 71 |
| 16 | 51 | 73 | 61 | B | 33 | 59 | 44/55 | 72 |
| 18 | 52 | 82 | 62/81 | C | 35 | 66 | 54 | |
| 26 | 53 | 6 | 67 | U | 39 | 68 | 61 | 84 |
| 31 | 56 | 11 | 69 | 42 | 45 | 73 | 62/81 | |
| | B | 40 | 70 | 43 | 51 | 82 | 67 | |



Genotypes

High-risk: 16, 18, 26, 31, 33, 35, 39, 45, 51, 52, 53, 56, 58, 59, 66, 68, 73, 82.

Low-risk: 6, 11, 40, 42, 43, 44, 54, 55, 61, 62, 67, 69, 70, 71, 72, 81, 84.



Sample types:

- ✓ Liquid-based cytology
- ✓ Paraffin-embedded tissues
- ✓ Cytology swabs



Single lyophilised master mix ready for use



Uses a direct sample: DNA extraction is not necessary

REFERENCE:
MAD-003930MU-HS24-24
MAD-003930MU-HS24-48

FORMAT:
24 TEST
48 TEST



REFERENCE: MAD-003938M-HS | 24 TEST

STD DIRECT FLOW CHIP KIT

Simultaneous detection of 12 pathogens that cause sexually transmitted diseases



Pathogens

Chlamydia trachomatis (serovars A-K),

Chlamydia trachomatis (serovars L1-L3),

Haemophilus ducreyi,

Herpes simplex virus 1,

Herpes simplex virus 2,

Mycoplasma genitalium,

Mycoplasma hominis,

Neisseria gonorrhoeae,

Treponema pallidum,

Trichomonas vaginalis,

Ureaplasma (*urealyticum/parvum*).



Sample types:

- ✓ Urine.
- ✓ Semen.
- ✓ Endocervical, pharyngeal, perianal, urethral and vaginal swabs.
- ✓ Liquid-based cytology.




Single lyophilised master mix ready for use.



Uses a direct sample: DNA extraction is not necessary.

| | | | | | | | | |
|----|-------------|-------|------|----|-------------|-------|------|--|
| B | MG | | | | CTS-1 | | B | |
| B | | MH | | | | CTS-2 | | |
| CI | TV | | UU-P | | MG | | UU-P | |
| BG | | HD | | | | MH | | |
| | HSV-1/HSV-2 | | NG | B | TV | | NG | |
| | | HSV-1 | | CI | | HD | | |
| | CTS-1 | | TP | BG | HSV-1/HSV-2 | | TP | |
| | | CTS-2 | | | | HSV-1 | | |
| | B | | pcCT | | | | pcCT | |



REFERENCE: MAD-003946M-HS | 24 TEST

MDR FLOW CHIP KIT

Simultaneous detection
of 5 bacterial species
and 56 resistance markers



Sample types:

- ✓ Blood cultures.
- ✓ Rectal and nasal swabs and rectal and nasal swabs as a single sample.
- ✓ Bacterial colonies.



Two lyophilised
master mixes
ready for use.



Uses a direct
sample: DNA
extraction is not
necessary.

| Resistance gene | Subclass |
|-----------------------------------------------------|-------------------------------------|
| mecA | BETA-LACTAM |
| mecC | BETA-LACTAM |
| VanA | VANCOMYCIN |
| VanB | VANCOMYCIN |
| blaSHV* (detección de mutaciones fenotipo BLEE) (3) | CEPHALOSPORIN |
| blaCTX* | CEPHALOSPORIN |
| KPC* | CARBAPENEM |
| SME* | CARBAPENEM |
| NMC/IMI* | CARBAPENEM |
| GES* | CEPHALOSPORIN/CARBAPENEM |
| VIM* | CARBAPENEM |
| OXA_23_like* | CARBAPENEM |
| SPM* | CARBAPENEM |
| OXA_51_like* | CARBAPENEM |
| NDM* | CARBAPENEM |
| OXA_58_like* | CARBAPENEM |
| OXA_24_like* | CARBAPENEM |
| GIM* | CARBAPENEM |
| OXA_48_like* | CARBAPENEM |
| SIM* | BETA-LACTAM |
| IMP_3, 15, 19_like* | CARBAPENEM |
| ctr | MACROLIDE/LINCOSAMIDE/STREPTOGRAMIN |
| mcrA | MACROLIDE |
| mef | MACROLIDE |
| ermA | MACROLIDE |
| ermB | MACROLIDE |
| ermC | MACROLIDE |
| aac(6)-Ib | AMIKACIN/KANAMYCIN/TOBRAMYCIN |
| armA | GENTAMICIN |

| Resistance gene | Subclass |
|--------------------------------------------------------------------------|--------------------|
| rmtB* | AMINOGLYCOSIDE |
| rmtC | AMINOGLYCOSIDE |
| rmtF* | AMINOGLYCOSIDE |
| DHA* | CEPHALOSPORIN |
| CMY* | CEPHALOSPORIN |
| sul1 | SULFONAMIDE |
| sul2 | SULFONAMIDE |
| sul3 | SULFONAMIDE |
| Point mutations in the <i>E. coli</i> DNA gyrase A subunit (4) | QUINOLONE |
| Point mutations in the <i>P. aeruginosa</i> DNA gyrase A subunit (3) | |
| Point mutations in the <i>E. coli</i> DNA topoisomerase IV A subunit (2) | |
| qnrA* | QUINOLONE |
| qnrB* | QUINOLONE |
| qnrS* | QUINOLONE |
| oqxA* | PHENICOL/QUINOLONE |
| oqxB* | PHENICOL/QUINOLONE |
| mcr1* | COLISTIN |
| mcr2* | COLISTIN |
| catB3 | CHLORAMPHENICOL |

Total: 56 resistance markers. *generic detection of these markers covering most allelic variants described to date

Bacteria

Staphylococcus aureus
Escherichia coli
Klebsiella pneumoniae
Pseudomonas aeruginosa
Acinetobacter baumannii

| | | | | | | | |
|----|-------|------|-----|----|-------|------|-----|
| B | | | | | | | B |
| B | | | | | | | MTB |
| CI | NEISS | AGAL | TPA | | LIS | CRYP | BOR |
| BG | | | | | | | |
| | SPNEU | HINF | COX | B | NEISS | AGAL | TPA |
| | | | | CI | | | |
| | LIS | CRYP | BOR | BG | SPNEU | HINF | COX |
| | MTB | | | | | | |
| | B | | | | | | |

REFERENCE: MAD-003935M-HS | 24 TEST

BACTERIAL CNS FLOW CHIP KIT

Simultaneous detection of
bacteria and fungi that cause
central nervous system infections



Bacteria and fungi

Mycobacterium tuberculosis complex, *Streptococcus pneumoniae*, *Streptococcus agalactiae*, *Haemophilus influenzae*, *Listeria monocytogenes*, *Treponema*

pallidum, *Neisseria meningitidis*, *Coxiella burnetii*, *Borrelia burgdorferi* and *Cryptococcus neoformans*.




Sample types:

- ✓ Purified DNA from cerebrospinal fluid.



Single master mix
ready for use.



REFERENCE: MAD-003936M-HS | 24 TEST SEPSIS FLOW CHIP KIT

Simultaneous detection of
36 bacterial species and
20 resistance markers

Gram+ bacteria

Streptococcus spp.
Streptococcus pneumoniae
Streptococcus pyogenes
Streptococcus agalactiae
Coagulase-negative staphylococcus
Staphylococcus aureus
Listeria monocytogenes
Enterococcus spp.
Neisseria meningitidis

Fungus

Candida spp.
Candida albicans

Gram- bacteria

Acinetobacter baumannii
Serratia marcescens
Klebsiella pneumoniae
Escherichia coli
Pseudomonas aeruginosa
Enterobacter spp.
Proteus spp./ *Morganella morganii*

Resistance

Methicillin: mecA
Vancomycin: vanA and vanB
Carbapenemase: KPC, SME, NMC/IMI
GES, VIM, GIM, SPM, NDM, SIM, IMP3, 15, 19_like, OXA23_like, OXA24_like, OXA48_like, OXA51_like, OXA58_like
β-lactamase: SHV and CTX-M



Sample types:

- ✓ Blood cultures.
- ✓ Rectal swabs.
- ✓ Bacterial colonies.



Two lyophilised master mixes ready for use.



Uses a direct sample: DNA extraction is not necessary.

| | | | | | | | | | |
|--------|-----------|----------|---------|-------|--------|-----------|----------|---------|-------|
| B | | LIS | kpc | spm | | ECOLI | vanB | | B |
| B | ABAU | ENTEROC | sme | ndm | | ENTEROB | vanA | ges | oxa23 |
| CI | SMAR/KLEB | PAER | nmc/imi | sim | | | mecA | vim | oxa24 |
| BG | SAGAL | KLEB | SPYOG | imp | SMALTO | CALB | | gim | oxa48 |
| | STAPHYL | STREP | blaSHV | | CAND | | PROT/MOR | kpc | oxa51 |
| SPENU | SA | NEISS | blaCTX | | B | ABAU | LIS | spm | oxa58 |
| | ECOLI | PROT/MOR | ges | oxa23 | CI | SMAR/KLEB | ENTEROC | sme | ndm |
| SMALTO | ENTEROB | | vim | oxa24 | BG | SAGAL | PAER | mmc/imi | sim |
| CAND | | mecA | gim | oxa48 | | STAPHYL | KLEB | SPYOG | imp |
| | CALB | vanA | | oxa51 | SPNEU | SA | STREP | blaSHV | |
| | B | vanB | | oxa58 | | | NEISS | blaCTX | |



REFERENCE: MAD-003937M-HS | 24 TEST AMR DIRECT FLOW CHIP KIT

Detection of 20 antibiotic
resistance genes present
in Gram-negative and
Gram-positive bacteria



Sample types:

- ✓ Blood cultures.
- ✓ Rectal/nasal swabs and rectal and nasal swabs as a single sample.
- ✓ Bacterial colonies.



Single lyophilised master mix ready for use.



Uses a direct sample: DNA extraction is not necessary.

Organism/Resistance

Staphylococcus aureus
Methicillin: mecA
Vancomycin: vanA and vanB
Carbapenemase: KPC, SME, NMC/IMI
GES, VIM, GIM, SPM, NDM, SIM, IMP3, 15, 19_like, OXA23_like, OXA24_like, OXA48_like, OXA51_like, OXA58_like
β-lactamase: SHV and CTX-M

| | | | | | | | | | |
|----|----|------|---------|------------|----|----|------|-----------|------------|
| B | | | kpc | spm | | | vanB | blaSHV-S | B |
| B | | | sme | ndm | | | vanA | ges | oxa23_like |
| CI | | | nmc/imi | sim | | | mecA | vim | oxa24_like |
| BG | | | | imp_like | | | | gim | oxa48_like |
| | | | blaSHV | blaSHV-S | | | | kpc | oxa51_like |
| | SA | | blaCTX | blaSHV-SK | B | | | spm | oxa58_like |
| | | | ges | oxa23_like | CI | | | sme | ndm |
| | | | vim | oxa24_like | BG | | | nmc/imi | sim |
| | | mecA | gim | oxa48_like | | | | blaSHV-SK | imp_like |
| | | vanA | | oxa51_like | | SA | | blaSHV | |
| | B | vanB | | oxa58_like | | | | blaCTX | |

REFERENCE: MAD-003939M-HS | 24 TEST

RESPIRATORY FLOW CHIP KIT

Simultaneous detection of the main agents
(13 pathogens) that cause acute respiratory infections



Pathogens

Influenza virus: FluA, FluA H1N1 2009, FluA H3 y FluB, Adenovirus, Bocavirus, Coronavirus (including SARS-CoV-2), Metapneumovirus, Parainfluenza virus (types 1, 2, 3 and 4),

Respiratory syncytial virus (type A and type B), Rinovirus, Enterovirus, *Bordetella pertussis*, *Bordetella parapertussis*, *Mycoplasma pneumoniae*.




Sample types:

✓ Purified RNA/DNA from bronchial lavages and nasopharyngeal swabs.



Two lyophilised master mixes ready for use.

| | | | | | | | | |
|--------|-----------|----------|----------|----------|-----------|-------|----------|--|
| B | FluA | PIV-1 | CoV-OC43 | | RNaseP | RSV-A | B | |
| B | FluA-H1N1 | PIV-2 | BP | | BG | RSV-B | CoV-229E | |
| CI-1 | FluA-H3 | PIV-3 | BPP | | | RhV | CoV-HKU1 | |
| CI-2 | FluB | PIV-4 | MP | | | PIV-1 | CoV-NL63 | |
| RNaseP | MPV | AdV | EV | B | FluA | PIV-2 | BPP | |
| BG | RSV-A | Bov | CoV-2 | CI-1 | FluA-H1N1 | PIV-3 | MP | |
| | RSV-B | CoV-229E | SARS | CI-2 | FluA-H3 | PIV-4 | EV | |
| | RhV | CoV-HKU1 | | CoV-OC43 | FluB | AdV | CoV-2 | |
| | B | CoV-NL63 | | BP | MPV | Bov | SARS | |



REFERENCE: MAD-003940M-HS | 24 TEST

TICK-BORNE BACTERIA FLOW CHIP KIT

Simultaneous detection of 7
arthropod-borne bacterial species



Bacteria

Rickettsia,
Rickettsia Typhus group,
Rickettsia Spotted fever group,
Erlichia,

Anaplasma,
Francisella,
Bartonella,
Borrelia,
Coxiella.



Sample types:

✓ Purified RNA/DNA from serum, biopsies, blood, cerebrospinal fluid, arthropods that carry the bacteria in question and animal reservoir samples.



Two lyophilised master mixes ready for use.

| | | | | | | | | |
|----|-----|--|-------|----|-----|--|-----|-------|
| B | | | FR | | | | B | |
| B | GR | | | | EH | | BOR | |
| CI | TG | | BAR | | AN | | | |
| BG | SFG | | | | | | COX | |
| | | | BAR-2 | B | | | | FR |
| | | | | CI | | | | |
| | EH | | BOR | BG | GR | | | BAR |
| | AN | | | | TG | | | BAR-2 |
| | B | | COX | | SFG | | | |



C/ Luis Fuentes Bejarano 60 Ed. Nudo Norte
(local 3) 41020 Seville - Spain

Phone: +34 954 93 32 00
pedidos@vitro.bio
www.vitro.bio