



FLU-COVID RT-PCR


REF: MAD-003942M

Kit for the detection of Influenza A (FluA),
Influenza B (FluB) and Severe Acute
Respiratory Syndrome-related Coronavirus 2
(SARS-CoV-2) by Multiplex One-Step Real-Time RT-PCR

 **vitro**
master diagnóstica®

CE IVD

FEATURES



Fast results



Target genes: SARS-CoV-2 specific N gene (CDC N2 region), FluA Matrix Protein and FluB Matrix Protein



**A single Ready to Use Master Mix
100 test
Lyophilized format available**



From extracted RNA
✓ Naso and oropharyngeal swabs
✓ Bronchoalveolar lavage (BAL)



Internal control and positive control included



Specificity and sensitivity of 100%



Characterized according to WHO and CDC recommendations



Limit of detection of 10 copies/reaction for all the three pathogens.



UNG enzyme incorporated to avoid contamination



Open system compatible with the most common extraction systems and RT-PCR platforms of the market